



# Robertson Microlit Laboratories

1705 US Highway 46/ Suite 1D / Ledgewood, N.J. 07852  
(973)966-6668 / Fax: (973)966-0136  
results@robertson-microlit.com www.robertson-microlit.com

## THERMAL ANALYSIS REQUEST

|   |  |
|---|--|
| SAMPLE #:   | CHARGE TO:   |
| DATE:   |  |
| REPORT TO:  | PO# OR CREDIT CARD INFO:   |
| COMPANY NAME:   | PHONE: ( ) -   |
| ADDRESS:  | FAX: ( ) -   |
| EMAIL: _____  |  |
| <b>ALL contact information mandatory</b>  |  |
| <b>DSC:</b> (10 mg sample required) <input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> cGMP (additional surcharge) |  |
| Use Default Parameters: <input type="checkbox"/>  | Use Client Parameters: <input type="checkbox"/>  |
| Sealed pan  | Type of pan:    sealed <input type="checkbox"/> perforated <input type="checkbox"/>  |
| Heating rate 10 <sup>0</sup> C/min  | Heating rate: _____  |
| Melting Point _____ <sup>0</sup> C  | Temp Range: _____ <sup>0</sup> C to _____ <sup>0</sup> C   |
| (temp range will be based on MP point)  | (normal range within 0 <sup>0</sup> C to 600 <sup>0</sup> C)   |
| ***If no MP provided, we will run melting point for DSC purposes (\$50 charge)  |  |
| Type of transition expected (if known): _____   |  |
| <b>TGA:</b> (10 mg sample required) <input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> cGMP (additional surcharge) |  |
| Use Default Parameters: <input type="checkbox"/>  | Use Client Parameters: <input type="checkbox"/>  |
| Heating rate 10 <sup>0</sup> C/min  | Heating rate: _____  |
| Temp range: ambient - 350 <sup>0</sup> C  | Temp Range: _____ <sup>0</sup> C to _____ <sup>0</sup> C   |
| (normal range within ambient to 700 <sup>0</sup> C)   |  |
| Graphic report options: <input type="checkbox"/> %weight loss vs. temp.   |  |
| <input type="checkbox"/> other _____  |  |
| <b>MELTING POINT:</b> <input type="checkbox"/> YES <input type="checkbox"/> NO    Amount of Sample Provided: _____                                |  |
| <input type="checkbox"/> cGMP (additional surcharge)  |  |
| <b>CUSTOMER SPECIAL INSTRUCTIONS</b>  | <b>LIST KNOWN HAZARDS and/or SAMPLE HANDLING PRECAUTIONS</b><br>(Samples WILL NOT be accepted if this section is left blank) |